

# Comprehensive liquid biopsy profiling enabled by PGDx elio plasma™ complete to facilitate precision oncology through decentralized access to testing

Kenneth C. Valkenburg<sup>1</sup>, Vito Caropreso<sup>1</sup>, Jesse Fox<sup>1</sup>, Christopher Gault<sup>1</sup>, Andrew Georgiadis<sup>1</sup>, Kelly M. R. Gerding<sup>1</sup>, Aanavi Karandikar<sup>1</sup>, Cynthia Maddox<sup>1</sup>, Paul McGregor<sup>1</sup>, David Riley<sup>1</sup>, James White<sup>1</sup>, Vuna Fa<sup>1</sup>

<sup>1</sup>Personal Genome Diagnostics, Inc., 3600 Boston Street, Baltimore, MD 21224, USA

## ABSTRACT


### Background

Despite the growing body of clinical evidence to support the utility of comprehensive genomic profiling (CGP) for advanced cancer patients, only a small fraction of individuals receive precision oncology guided treatment strategies. We set out to develop the PGDx elio™ plasma complete assay in order to guide both translational biomarker discovery, monitoring response, and CGP-informed precision oncology strategies in cell free DNA (cfDNA) from the plasma of cancer patients.

### Specific Objectives

- To develop a 500+ gene hybrid-capture next generation sequencing (NGS) assay as a decentralized kit for cell-free DNA from liquid biopsies
- To enable comprehensive detection of SNVs, indels, amplifications, translocations, MSI, and blood TMB
- To identify clinically relevant cancer variants, hotspots, and genomic signatures in plasma cfDNA
- To build an assay with high specificity, sensitivity, accuracy, precision, reproducibility, and robustness
- To enable compatibility on the NovaSeq 6000 Sequencer
- To build manual and automated versions of the assay for user flexibility

## ASSAY WORKFLOW



### General Assay Specifications for PGDx elio plasma complete

Parameter	Specification
<b>Panel Size</b>	2.1 MB
<b>Panel Content and Variant Type</b>	521 genes for SNV & Indels 38 genes for amplifications 21 genes for translocations MSI Status bTMB (Muts/Mb) LOH Status
<b>Reportable Range</b>	SNVs & indels: ≥ 0.1% VAF Translocations: ≥ 2 fusion reads Amplifications: ≥ 1.15-fold
<b>Sample Requirements</b>	Plasma cfDNA
<b>DNA Input Requirements</b>	25 ng recommended; 10 ng minimum
<b>Sample Pass Rate</b>	97.8% overall pass rate (271/277)
<b>Sequencing Platform/ Flow Cell</b>	NovaSeq 6000/S2 Flow cell
<b>Sequence Run</b>	2 x 150 bp
<b>Cases per Sequencing Run</b>	16
<b>Workflow</b>	Manual and Automated Available
<b>Variant Reporting</b>	Automated analytical pipeline and variant report generation
<b>Average Total Coverage</b>	~20,000x
<b>Average De-duplicated Error-corrected Coverage</b>	~2300x

## METHODS

- Specificity was assessed using plasma cfDNA from 20 noncancerous donors at 25 ng input.
- Sensitivity was assessed in 3 series of dilutions of variant-positive DNA into wild-type DNA to multiple levels above and below the theoretical limit of detection
- Accuracy was assessed by comparing reported variant calls in clinical plasma samples and reference materials to orthogonal assays in 71 samples (including bladder, colorectal, gastric, lung, ovarian, and thyroid cancer) at 25 ng input.
- Precision, reproducibility, and repeatability was assessed in 2 well-characterized samples at 25 ng input in triplicate across 2 operators and 3 non-consecutive runs.
- DNA input was assessed at multiple inputs across 350 clinical samples. In addition, concordance between the 25 ng recommended input and the 10 ng minimum input was assessed in 15 samples.
- Comparison of manual and automated lab assays: 45 samples were processed manually and on the Tecan Freedom EVO 150.

## RESULTS

### Analytical Accuracy Primary Endpoint Results Compared to Targeted NGS Panels (n=64)

Analyte	PPA	NPA
SNVs	92.7%	99.9%
Indels	94.4%	99.9%
Translocations	82.4%	100%
Amplifications	89.3%	96.4%
MSI	100%	100%
bTMB	0.72 Spearman Correlation Coefficient	


### Analytical Accuracy Primary Endpoint Results Compared to Competitor 500+ Gene cfDNA Assay (n=7)

Analyte	PPA (2-sided 95% CI)	NPA (2-sided 95% CI)
SNVs	92.2% (95/103)	99.99% (8414930/8414933)
Indels	83.3% (10/12)	99.99% (8415023/8415024)
Translocations*	50% (2/4)	99.7% (366/367)
Amplifications†	76% (19/25)	98.6% (141/143)
MSI	N/A	100% (7/7)
bTMB	0.73 Spearman Correlation Coefficient	

\* Discordant calls: BRAF-ZC3HAV1 non-actionable, low fusion read count; EWSR1-RP11-9L18.2, filtered due to fusion with pseudogene  
† All calls were < 2-fold except 1 concordant MET call

### SNV and Indel LoD Summary

Variant Category	Number of Variants	Observed Range <sup>2</sup>	Median LoD <sup>3</sup>
Clinically Relevant SNVs & Indels	10 (9 SNV; 1 indel)	0.32% - 0.78%	0.40% VAF
Panel-wide SNVs & Indels	263 (245 SNV; 18 indels)	0.34% - 1.75%	1.16% VAF
Translocations	2	0.22% - 0.82%	0.24% and 0.41% FRF
Amplifications	1	1.28 - 1.30-fold	1.29-fold



### Analytical Specificity: Limit of Blank Primary Endpoint Results

Analyte Assessed	Specificity (n/N)
SNVs (clinically relevant)	100% (4260/4260)
Panel-Wide SNVs	99.9999% (28009535/28009540)
Indels (clinically relevant)	100% (1780/1780)
Panel-wide Indels	99.9999% (28009528/28009540)
Translocations	100% (420/420)
Amplifications	100% (760/760)
MSI	100% (20/20)
bTMB*	100% (20/20)


\*Confirmation that non-cancerous samples bTMB score was below the established Limit of Blank of 1.0 Muts/Mb

## RESULTS


### Manual vs. Automated Assay Workflow

Variable Assessed	Manual	Automated
<b>Library Preparation Yield</b>	All samples had ≥ 108 ng/μL yield	
<b>Total Coverage</b>	Mean = 7.3% CV	
<b>De-duplicated error-corrected coverage</b>	Mean = 5.9% CV	
<b>GC Bias</b>	Mean = 0.42% CV	
<b>Pass/Fail Status</b>	100% success rate	
<b>Clinically Relevant SNV &amp; Indel Concordance</b>	98.4% PPA	99.94% NPA
<b>Panel-wide SNV &amp; Indel Concordance</b>	98.0% PPA (≥1.0% VAF)	92.2% PPA (all variants assessed)
<b>Amplification Concordance</b>	95.79% PPA	99.94% NPA
<b>Translocation Concordance</b>	100% PPA	99.78% NPA
<b>MSI Concordance</b>	100% PPA	100% NPA
<b>bTMB Concordance</b>	Mean = 6.5% CV	


### MSI Reproducibility



### bTMB Reproducibility



### De-duplicated Error-corrected Coverage per DNA Input



DNA Input Bin (ng)	N	Observed Pass Rate
10 - 19.9	68	85.3%
20-24.9	51	92.9%
25	231	98.2%
Total	350	

### Precision and Reproducibility: Average Positive Agreement (APA) and Average Negative Agreement (ANA) Results

Alteration Type	Overall		Inter-Operator		Inter-Day		Inter-Run		Within-Run	
	APA (%)	ANA (%)	APA (%)	ANA (%)	APA (%)	ANA (%)	APA (%)	ANA (%)	APA (%)	ANA (%)
<b>Alteration Type</b>	95% CI	95% CI	95% CI	95% CI	95% CI	95% CI	95% CI	95% CI	95% CI	95% CI
<b>MSI</b>	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%
<b>Amplifications (FGF4, FGF19, CCND1)</b>	95.9%	95.9%	93.4%	93.4%	90.4%	90.4%	90.4%	90.4%	90.4%	90.4%
<b>Translocations (ALK, RET)</b>	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%	100.0%
<b>Clinically Relevant SNVs and Indels</b>	99.5%	99.6%	99.1%	99.3%	98.7%	98.9%	98.7%	98.9%	98.7%	98.9%
<b>Panel-wide SNVs and Indels</b>	92.5%	94.0%	92.5%	93.9%	92.7%	94.1%	92.7%	94.1%	92.7%	94.1%
<b>Panel-wide SNVs and Indels</b>	93.5%	94.8%	93.7%	94.9%	94.1%	95.3%	94.1%	95.3%	94.1%	95.3%

## CONCLUSIONS

PGDx built a 500+ gene decentralized NGS panel for plasma with high sensitivity, specificity, and reproducibility with manual and automated lab assay options. Through pre-defined, robust design control processes, the product was designed, verified, and validated. The PGDx elio plasma complete assay was designed for ease of use both in the manual setting and automated setting (using the Tecan Freedom EVO 150), with a multiplex configuration that reduces lab steps and time. The assay is run with 16 samples per kit from library preparation through sequencing over the course of ≤ 3.5 days. While the recommended DNA input into the assay is 25 ng, the assay performs well down to 10 ng input. The NovaSeq instrument provides sufficient sequencing depth of coverage for high sensitivity over a large panel, which includes a genomic region of 2.1 Mb covering 521 genes. The automated software produces in-depth reports after an average of 3.5 days from completion of sequencing. The PGDx elio plasma complete reports provide information on single nucleotide variants (SNVs) and insertions and deletions (indels) for all 521 genes on the panel, as well as being able to report amplifications in 38 genes and translocations in 21 genes (see Appendix). The product also reports microsatellite instability (MSI) status, loss of heterozygosity (LOH) status, and a numerical value for blood tumor mutation burden (bTMB) for every case. Overall, elio plasma complete is a highly sensitive, specific, and reproducible plasma assay.